

WHO-GMP, ISO 9001:2015, ISO 13485:2016 CERTIFIED COMPANY


We Perform to Lead



Division of - Gajham Exports Pvt. Ltd.
Government Recognised Export House

FULLY MOTORISED BED

EXCELLENZA



Overall Dimension : 2230mm (L) x 1050mm (W)
Suitable for mattress : 1980mm (L) x 880mm (W)

Minimum Height : 450mm (without mattress)
Maximum Height : 750mm (without mattress)

Features:

- 60mm x 30mm CRCA tubular frame work.
- Four section perforated CRCA sheet top.
- Height adjustment, Backrest, Kneerest, Trendelenburg & R T, Electrically Operated by Handset.
- CPR Function, MCP Box & Battery Back-up.
- Removable Polymer Head and Foot Boards.
- Split type Polymer Railing (set of 4).
- S.S. Telescopic Saline Rod with two locations.
- Central & directional locking castors.
- Four side corner buffers.
- Urine bag holder.
- Pretreated and Epoxy powder coated.

• OPTIONAL ACCESSORIES AT EXTRA COST:

- a) Infusion pump rod with telescopic I.V. Rod.
- b) PU foam mattress.
- c) Bed extension facility.



47 (B) SEMI MOTORISED BED



Overall Dimension : 2230mm (L) x 1050mm (W)
Suitable for mattress : 1980mm (L) x 880mm (W)

Minimum Height : 450mm (without mattress)
Maximum Height : 750mm (without mattress)

Features:

- 60mm x 30mm CRCA tubular frame work.
- Four section perforated CRCA sheet top.
- Height adjustment, Backrest & Knee rest is Electrically Operated with Handset.
- Removable Polymer Head and Foot Board.
- Split type Polymer Railings (set of 4).
- S.S. Telescopic Saline Rod with two locations.
- Four heavy duty swivel castors of 125mm Dia (two with brake).
- Four side corner buffers.
- Pretreated and Epoxy powder coated.

• OPTIONAL ACCESSORIES AT EXTRA COST:

- Nursing Control Panel.
- CPR function.
- Trendelenburg & Reverse Trendelenburg.
- Central & directional locking castors.
- Infusion pump Rod with telescopic I.V. Rod.
- PU foam mattress.
- Urine bag holder.
- Bed extension facility.



ACUTE CARE (DELUXE)

- 08 (E) ACUTE CARE DELUXE
HEIGHT ADJUSTABLE ICU BED WITH SPLIT TYPE SIDE RAILINGS



Overall Size : 2350mm (L) x 1050mm (W) x Ht. adj. from 550mm to 750mm.
Table size for mattress : 1980mm (L) x 890mm (W)

Features:

- 60mm x 30mm CRCA tubular frame work.
- Four section perforated CRCA sheet top.
- Height adjustment, Backrest, Kneerest, Trendelenburg & Reverse Trendelenburg operated by screw and lever mechanism with four separate attached S.S. folding handles.
- Polymer moulded Head and Foot Board.
- Split type side Railings (set of 4).
- S.S. Telescopic Saline Rod with four locations.
- Four heavy duty swivel castors of 125mm Dia. (two with brake).
- Four side corner buffers.
- Pretreated and Epoxy powder coated.

• OPTIONAL ACCESSORIES AT EXTRA COST:

- a) Central & Directional Locking Castors
- b) Infusion pump Rod with telescopic I.V. Rod.
- c) Four section PU foam Mattress.
- d) Urine Bag Holder.



ACUTE CARE (DELUXE)

- 08 (F) ACUTE CARE DELUXE
HEIGHT ADJUSTABLE ICU BED WITH POLYMER RAILINGS (SET OF TWO)



Overall Size : 2350mm (L) x 1050mm (W) x Ht. adj. from 550mm to 750mm.
Suitable size for mattress : 1980mm (L) x 890mm (W)

Features:

- 60mm x 30mm CRCA tubular frame work.
- Four section perforated CRCA sheet top.
- Height adjustment, Backrest, Kneerest, Trendelenburg & Reverse Trendelenburg operated by screw and lever mechanism with four separate attached S.S. folding handles.
- Polymer moulded Head and Foot Boards.
- A pair of polymer moulded side railings.
- S.S. Telescopic Saline Rod with two locations.
- Four heavy duty swivel castors of 125mm Dia. (two with brake).
- Four side corner buffers.
- Pretreated and Epoxy powder coated.

• OPTIONAL ACCESSORIES AT EXTRA COST:

- a) Central & Directional Locking Castors.
- b) Infusion pump Rod with telescopic I.V. Rod.
- c) Four section PU foam Mattress.
- d) Urine bag holder.



ACUTE CARE (DELUXE)

08 (D) ACUTE CARE DELUXE HEIGHT ADJUSTABLE ICU BED WITH NEW DESIGN COLLAPSIBLE SIDE RAILINGS



Size : 2350mm (L) x 1050mm (W) x Ht. adj. from 550mm to 750mm.
Suitable size for mattress : 1980mm (L) x 890mm (W)

Features:

- 60mm x 30mm CRCA tubular frame work.
- Four section perforated CRCA sheet top. Backrest, Kneerest, Trendelenburg, Reverse Trendelenburg & Height adjustment operated by screw and lever mechanism with four separate attached S.S. folding handles.
- Polymer moulded Head and Foot Board.
- A pair of new design collapsible railings (1320mm length).
- S.S. Telescopic Saline Rod with four locations.
- Four heavy duty swivel castors of 125mm Dia. (two with brake).
- Four side corner buffers.
- Pretreated and Epoxy powder coated.


• OPTIONAL ACCESSORIES AT EXTRA COST:

- a) Central & Directional Locking Castors.
- b) Infusion pump rod with telescopic I.V. Rod.
- c) Four section PU foam Mattress.
- d) Urine bag holder.



ACUTE CARE (DELUXE)

07 (D) ACUTE CARE DELUXE FIXED HEIGHT ICU BED WITH NEW DESIGN COLLAPSIBLE SIDE RAILINGS



Overall Size : 2350mm (L) x 1050mm (W) x 590mm (H).
Suitable size for mattress : 1980mm (L) x 890mm (W)

Features:

- 60mm x 30mm CRCA tubular frame work.
- Four section perforated CRCA sheet top.
- Backrest, Kneerest, Trendelenburg & Reverse Trendelenburg operated by screw and lever mechanism with three separate attached S.S.folding handles.
- Polymer Head and Foot Boards.
- A pair of new design collapsible side railings (1320mm length).
- S.S. Telescopic Saline Rod with four locations.
- Four heavy duty swivel castors of 125mm Dia. (two with brake).
- Four side corner buffers.
- Pretreated and Epoxy powder coated.


• OPTIONAL ACCESSORIES AT EXTRA COST:

- a) Central & Directional Locking Castors.
- b) Infusion pump Rod with telescopic I.V. Rod.
- c) Four Section PU foam mattress.
- d) Urine bag holder.



ACUTE CARE (RECOVERY)

05 (A) ACUTE CARE RECOVERY BED
WITHOUT RAILINGS & CASTORS



05 (B) ACUTE CARE RECOVERY BED
WITH SWING TYPE SIDE RAILINGS & CASTORS



Overall Size 1) : 2250mm (L) x 1020mm (W) x 590mm (H).
Suitable size for mattress : 1905mm (L) x 890mm (W)

Overall Size 2) : 2170mm (L) x 930mm (W) x 590mm (H).
Suitable size for mattress : 1830mm (L) x 810mm (W)

Features:

- 60mm x 30mm CRCA tubular frame work.
- Two section perforated CRCA sheet top.
- Backrest, Trendelenberg and Reverse Trendelenberg operated by screw and lever mechanism with attached S.S. folding handles.
- 32mm Dia. Stainless Steel bows with colored metal panel.
- A pair of swing away side railings.
- Four locations for I.V. Rod.
- Four heavy duty swivel castors of 125mm dia (two with brake).
- Pretreated and Epoxy powder coated.

• OPTIONAL ACCESSORIES AT EXTRA COST:

- a) Saline rod.
- b) PU foam mattress.

ACUTE CARE

08 (B) ACUTE CARE HEIGHT ADJUSTABLE ICU BED WITH FULL LENGTH COLLAPSIBLE SIDE RAILINGS



Overall Size: 1) : 2350mm (L) x 1020mm (W) x Ht. adj. from 550mm to 750mm.
Suitable size for mattress : 1980mm (L) x 890mm (W)

Overall Size: 2) : 2250mm (L) x 930mm (W) x Ht. adj. from 550mm to 750mm.
Suitable size for mattress : 1905mm (L) x 810mm (W)

- OPTIONAL ACCESSORIES AT EXTRA COST:
- a) Infusion pump Rod with telescopic I.V. Rod.
- b) Four Section PU foam mattress.
- c) Urine bag holder.

Features:

- 60mm x 30mm CRCA tubular frame work.
- Four section perforated CRCA sheet top.
- Height adjustment, Backrest, Kneerest, Trendelenburg & Reverse Trendelenburg operated by screw and lever mechanism with four separate attached S.S. folding handles.
- 32mm Dia Stainless Steel bows with colored metal panel.
- A pair of Full Length collapsible side railings.
- S.S. Telescopic Saline Rod with four locations.
- Four heavy duty swivel castors of 125mm Dia. (two with brake).
- Four side corner Buffers.
- Pretreatment and Epoxy powder coated.

ACUTE CARE

08 (C) ACUTE CARE HEIGHT ADJUSTABLE ICU BED WITH NEW DESIGN COLLAPSIBLE SIDE RAILINGS



Overall Size: 1) : 2350mm (L) x 1020mm (W) x Ht. adj. from 550mm to 750mm.
Suitable size for mattress : 1980mm (L) x 890mm (W)

Overall Size: 2) : 2250mm (L) x 930mm (W) x Ht. adj. from 550mm to 750mm.
Suitable size for mattress : 1905mm (L) x 810mm (W)

Features:

- 60mm x 30mm CRCA tubular frame work.
- Four section perforated CRCA sheet top.
- Height adjustment, Backrest, Kneerest, Trendelenburg & Reverse Trendelenburg operated by screw and lever mechanism with four separate attached S.S. folding handles.
- 32mm Dia Stainless Steel bows with colored metal panel.
- New Design Collapsible Side Railings (1320 mm length).
- S.S. Telescopic Saline Rod with four locations.
- Four heavy duty swivel castors of 125mm Dia (two with brake).
- Four side Corner Buffers.
- Pretreated and Epoxy powder coated.

• OPTIONAL ACCESSORIES AT EXTRA COST:

- a) Infusion pump Rod with telescopic I.V. Rod.
- b) Four Section PU Foam Mattress.
- c) Urine bag holder.

ACUTE CARE

- 07 (A) ACUTE CARE BED FIXED HEIGHT WITH SWING AWAY RAILINGS



Overall Size 1) : 2350mm (L) x 1020mm (W) x 590mm (H)
Suitable size for mattress : 1985mm (L) x 890mm (W)
Overall Size 2) : 2250mm (L) x 930mm (W) x 590mm (H)
Suitable size for mattress : 1905mm (L) x 810mm (W)

Features:

- 60mm x 30mm CRCA tubular frame work.
- Four section perforated CRCA sheet top.
- Backrest, Kneerest, Trendelenberg and Reverse Trendelenberg operated by screw and lever mechanism with attached S.S. folding handles.
- 32mm Dia. Stainless Steel bows with colored metal panel.
- A pair of stainless steel swing away side railings.
- Four locations for I.V. Rod.
- Four heavy duty swivel castors of 125mm dia (two with brake)
- Four side corner buffers.
- Pretreated and Epoxy powder coated.

• OPTIONAL ACCESSORIES AT EXTRA COST:

- a) Infusion pump rod with telescopic I.V. Rod.
- b) PU foam mattress.
- c) Urine bag holder.

- 07 (B) ACUTE CARE BED FIXED HEIGHT WITH FULL LENGTH COLLAPSIBLE RAILINGS

Overall Size 1) : 2350mm (L) x 1020mm (W) x 590mm (H).
Suitable size for mattress : 1980mm (L) x 890mm (W)
Overall Size 2) : 2250mm (L) x 930mm (W) x 590mm (H).
Suitable size for mattress : 1905mm (L) x 810mm (W)

Features:

- 60mm x 30mm CRCA tubular frame work.
- Four section perforated CRCA sheet top.
- Backrest, Kneerest, Trendelenberg & Reverse Trendelenberg operated by screw and lever mechanism with three separate attached S.S. folding handles.
- 32mm Dia Stainless Steel bows with colored metal panel.
- S.S. Telescopic Saline Rod with four locations.
- A pair of Full Length Collapsible side railings
- Four heavy duty swivel castor of 125mm Dia (two with brake).
- Four Side Corner buffers.
- Pretreated and Epoxy powder coated.

• OPTIONAL ACCESSORIES AT EXTRA COST:

- a) Infusion pump Rod with telescopic I.V. Rod.
- b) PU foam mattress.
- c) Urine bag holder.



ACUTE CARE

- 07 (C) ACUTE CARE FIXED HEIGHT BED WITH NEW DESIGN COLLAPSIBLE SIDE RAILINGS



Overall Size 1) : 2350mm (L) x 1020mm (W) x 590mm (H).

Suitable size for mattress : 1980mm (L) x 890mm (W)

Overall Size 2) : 2250mm (L) x 930mm (W) x 590mm (H).

Suitable size for mattress : 1905mm (L) x 810mm (W)

Features:

- 60mm x 30mm CRCA tubular frame work.
- Four section perforated CRCA sheet top.
- Backrest, Kneerest, Trendelenberg & Reverse Trendelenberg operated by screw and lever mechanism with three separate attached S.S. folding handles.
- 32mm Dia Stainless Steel bows with colored metal panel.
- S.S. Telescopic Saline Rod with four locations.
- A pair of New Design Collapsible side railings (1320 mm length).
- Four heavy duty swivel castors of 125mm Dia (two with brake).
- Four side corner buffers.
- Pretreated and Epoxy powder coated.

• OPTIONAL ACCESSORIES AT EXTRA COST:

- Infusion pump Rod with telescopic I.V. Rod.
- PU foam mattress.
- Urine bag holder.



PRIMA (DELUXE)

04 (B) FOWLER BED DELUXE - PRIMA



Overall Size 1) : 2350mm (L) x 1050mm (W) x 590mm (H).
Suitable size for mattress : 1980mm (L) x 890mm (W)

Features:

- 60mm x 30mm CRCA tubular frame work.
- Four section perforated CRCA sheet top.
- Backrest and Kneerest operated by screw and lever mechanism with attached S.S. folding handle.
- Polymer moulded head and foot boards.
- A pair of new design collapsible side railings (1320mm length).
- Four heavy duty swivel castor of 125mm Dia. two with brake.
- Four side corner buffers.
- Four locations for I.V. Rod.
- Pretreated and Epoxy powder coated.

• OPTIONAL ACCESSORIES AT EXTRA COST:

- a) Saline Rod.
- b) PU foam mattress.

PRIMA (DELUXE)

03 (C) SEMI FOWLER BED (DELUXE) - PRIMA



Overall Size 1) : 2060mm (L) x 1050mm (W) x 590mm (H).
Suitable size for mattress : 1905mm (L) x 890mm (W)

Features:

- 60mm x 30mm CRCA tubular frame work.
- Two section perforated CRCA sheet top.
- Backrest operated by screw and lever mechanism with attached S.S. folding handle.
- Polymer moulded head and foot boards.
- A pair of new design collapsible side railings (1320mm length).
- Four heavy duty swivel castor of 125mm Dia. two with brake.
- Four side corner buffers.
- Four locations for I.V. Rod.
- Pretreated and epoxy powder coated.

• OPTIONAL ACCESSORIES AT EXTRA COST:

- a) Saline Rod.
- b) Urine bag holder.
- c) PU foam mattress.

PRIMA

04 (A) FOWLER BED (STANDARD)



Overall Size 1) : 2160mm (L) x 925mm (W) x 590mm (H)
Suitable size for mattress : 1980mm (L) x 890mm (W)

Overall Size 2) : 2060mm (L) x 850mm (W) x 590mm (H)
Suitable size for mattress : 1905mm (L) x 810mm (W)

Features:

- 60mm x 30mm CRCA tubular frame work.
- Four section perforated CRCA sheet top.
- Backrest and Kneerest operated by screw and lever mechanism with attached S.S. folding handles.
- 32mm Dia. Stainless Steel bows with colored metal panel.
- Four locations for I.V. Rod.
- Pretreated and Epoxy powder coated.

OPTIONAL ACCESSORIES AT EXTRA COST:

- a) Castors (set of four) two with brake.
- b) Saline Rod.
- c) A pair of side Railings.
- d) PU foam mattress.

03 (A) SEMI FOWLER BED (STANDARD)

Overall Size 1) : 2060mm (L) x 925mm (W) x 590mm (H).

Suitable size for mattress : 1905mm (L) x 890mm (W)

Overall Size 2) : 1985mm (L) x 850mm (W) x 590mm (H).

Suitable size for mattress : 1830mm (L) x 810mm (W)

Features:

- 60mm x 30mm CRCA tubular frame work.
- Two section perforated CRCA sheet top.
- Backrest operated by screw and lever mechanism with attached S.S. folding handle.
- 32mm Dia. Stainless Steel bows.
- Four locations for I V Rod.
- Pretreated and Epoxy powder coated.
- Sizes available as per users need.

OPTIONAL ACCESSORIES AT EXTRA COST:

- a) Castors (set of four) two with brake.
- b) Saline Rod.
- c) A pair of side railings.
- d) Crib attachment.
- e) Colored Metal Panel.
- f) PU foam mattress.
- g) Urine bag holder.

03 (B)



PRIMA

01 PLAIN HOSPITAL BED



Overall Size 1) : 2010mm (L) x 925mm (W) x 590mm (H).
Suitable size for mattress : 1905mm (L) x 890mm (W)

Overall Size 2) : 1935mm (L) x 850mm (W) x 590mm (H).
Suitable size for mattress : 1830mm (L) x 810mm (W)

Features:

- 60mm x 30mm CRCA tubular frame work.
- Perforated CRCA sheet top.
- 32mm dia Stainless Steel bows.
- Four locations for I.V.Rod.
- Pretreated and Epoxy powder coated.
- Sizes available as per users need.

• OPTIONAL ACCESSORIES AT EXTRA COST:

- a) Castors Set of four (two with brake).
- b) Saline Rod.
- c) A pair of side railings.
- d) Urine bag holder.
- e) Crib attachment.
- f) Colored metal bow panel.
- g) Mattress.

02 HOSPITAL BED WITH MANUAL BACKREST

Overall Size 1) : 2010mm (L) x 925mm (W) x 590mm (H).
Suitable size for mattress : 1905mm (L) x 890mm (W)

Overall Size 2) : 1935mm (L) x 850mm (W) x 590mm (H).
Suitable size for mattress : 1830mm (L) x 810mm (W)

Features:

- 60mm x 30mm CRCA tubular frame work.
- Two section perforated CRCA sheet top.
- Manually operated six position backrest on ratchet.
- 32mm dia Stainless Steel bows.
- Four locations for I.V.Rod.
- Pretreated and epoxy powder coated.
- Sizes available as per users need.

• OPTIONAL ACCESSORIES AT EXTRA COST:

- a) Castors Set of four (two with brake).
- b) Saline Rod
- c) A pair of side railings
- d) Crib attachment
- e) Colored metal bow panel.
- f) Urine bag holder.
- g) Mattress.




PRIMA

09 (B) ATTENDENT BED (SLIDING)

Overall Approx Size Open
1830mm (L) x 560mm (W) x 390mm (H)

Overall Approx Size Closed :
940mm (L) x 560mm (W) x 390mm (H)

- Sliding with folding mattress
- Pretreated and Epoxy powder coated



09 (A) STANDARD ATTENDENT

Overall Approx Size :
1830mm (L) x 610mm (W) x 470mm (H)
50mm x 25mm CRCA rectangular tubular frame work.

- Two locations for I.V. Rod.
- Pretreated and Epoxy powder coated,

06 PEADIATRIC BED WITH DROP SIDE RAILINGS

Overall size : L 1450mm x W 840 mm x H 590mm
Mattress Size : L1350mm x W 750mm

- Two locations for I.V. Rod
- Perforated CRCA Sheet Top
- Full length drop side Railings.
- Pretreated & Epoxy powder coated.



49 RECLINING CHAIR

Overall approx. Size : L 1730mm x w 558mm
Chair Position: L 780mm x W 558mm



WARD REQUISITE

12 BED SIDE LOCKER (STANDARD)

Overall approx Size : W 370 mm x D 305 mm x H 760 mm

Features:

- Two Stainless Steel tops with three side raised & front side down.
- Railings on three sides above main cabinet.
- Handle & magnetic latch on door.
- CRCA sheet box mounted on 20mm x 20mm CRCA tubes frames.
- Pretreated & Epoxy powder coated.



13 BED SIDE LOCKER (DELUXE)

Overall Approx Size : W 395 mm x D 370 mm x H 805 mm

Features:

- CRCA Sheet main cabinet and a drawer.
- Mounted on 50 mm Dia. castors.
- Stainless Steel top with three sides raised & front side down.
- Buffers at rear side.
- Pretreated & Epoxy powder coated.



14 BED SIDE STAND FULL S.S.

Overall Approx Size : W 390mm x D 360 mm x H 770 mm

Features:

- Three Stainless Steel Shelves.
- Mounted on Stainless Steel tubular frame.
- Optional : 50 mm Dia. Castors.

WARD REQUISITE

33 MONITOR TROLLEY

Overall size : W 965mm x D 660mm x H 1498mm
Shelves size : W 610mm x D 510mm

Features:

- Adjustable 5 shelves
- Mounted on 125mm Dia. Heavy duty castors 2 with brakes.
- Eight point electrical panel board with socket & switch.
- Co2 Cylinder cage.
- Available in Stainless Steel / Powder Coated.



A



B



A



B

38 INSTRUMENT CABINET

Overall Size : W 610mm x D 355mm x H 1535mm

Features:

- Rectangular tubular structure.
- Sheet Panel on back side.
- Five shelves, door & two side of 5mm thick fine glass.
- Available in Stainless Steel / Powder coated.

WARD REQUISITE



39 CRASH CART

Overall approx. Size : W 970mm x D 490mm x H 1630mm

Features:

- Six drawer modular system.
- Upper Two drawer with medicine containers.
- Detachable plastic bins.
- oxygen cylinder Cage .
- Top tray for monitor & pulse oxymeter.
- S.S. I.V. Rod & cardiac massage board.
- System mounted on 100 mm Dia castors two with brake.
- Available in Stainless Steel / Powder coated frame with S.S. top.

17 HEIGHT ADJUSTABLE BED SIDE TABLE

Top tray Size: L 760mm x W 360mm

Features:

- Height Adjustable from 740 mm to 1040 mm
- Stainless Steel top .
- Height adjustable by bevel gear mechanism .
- Mounted on 50mm Dia. castors.
- Pretreated & Epoxy powder coated.
- Optional : Height adjustment by gas spring.



MI - 45 WALL MOUNTING MONITOR STAND

Tray Size : W 270mm x D 380mm

- Pretreated & Epoxy powder coated.

O.T.EQUIPMENTS

23 (C) OPERATION TABLE (ALL FUNCTION)




Overall approx. size: L 1830 mm x W 710 mm x Ht. Adjustable from 735 mm to 970 mm

Features:

- Stainless steel top and stainless steel base.
- Kidney bridge.
- Removable & folding foot end with S.S. drainage tray.
- Height adjustment and Head low/up (T & RT) position by screw & lever mechanism.
- Lateral Tilt Position.
- Lithotomy crutches.
- Gear type angular arrangement for lithotomy bars.
- Mattress, Step Stool & Arm boards.
- Stainless Steel removable side railings

O.T.EQUIPMENTS



23 (A) HEIGHT AJUSTABLE OPERATION TABLE

Overall Size : 1830mm (L) x 710mm (W)

Features:

- Height adjustment from 735mm to 970mm.
- Stainless Steel top , Stainless Steel Base .
- Removable & folding foot end with S.S. drainage tray.
- Height adjustment and Head low/up (T & RT) position by screw & lever mechanism.
- Lithotomy bars with stirrups
- Mattress, Step Stool & Arm Boards.

• **OPTIONAL ACCESSORIES AT EXTRA COST:**

- Kidney bridge.
- Lithotomy crutches.

22 (A) SIMPLE OPERATION TABLE

Overall approx. size: L 1830 mm x W 710 mm x Ht 810 mm

Features:

- Stainless Steel top, Powder Coated Base.
- Removable & folding foot end with S.S. drainage tray.
- Head low/up (T & RT) position by screw & lever mechanism.
- Lithotomy Bars with stirrups & Arm Boards.
- Mattress & Step Stool.
- Also available in Complete Stainless Steel.



23 (D) OPHTHALMIC OPERATION TABLE

Features:

Overall approx size : L 1830 mm x W 690 mm x Hgt. adj. from 735mm to 970mm

- Stainless Steel Top, Stainless Steel Base.
- Adjustable head rest.
- Height adjustment & Head low/up (T & RT) position by screw & lever mechanism.
- Mattress & Step Stool.



GYNAEC CARE

27 TELESCOPIC LABOUR TABLE



Overall approx. Size : L 1830mm x W 870mm x H 800mm
Frame : 50mm x 25mm, Rectangular Stainless Steel Tube.

Features:

- Three Section Stainless Steel top, Stainless Steel Base.
- Backrest on ratchet.
- Hand grips on both side.
- Head low/up (T & RT) by screw mechanism.
- Middle section having perineal U cut with bowl.
- Sliding foot end.
- Railing on three sides.

• ACCESSORIES: Lithotomy bars, Stirrups, Step Stool & Mattress.

26 (B) LABOUR TABLE (FOOT END FOLDING)



26 (C) WITH BACKREST



Overall approx. Size : L 1830mm X W 870mm x H 800mm
Frame : 50mm x 25mm, Rectangular Stainless Steel Tube.

Features:

- Two Section Stainless Steel top & base with head low/up (T & RT) by screw mechanism
- Removable & folding foot end.
- Perineal 'U' cut with drainage tray.
- Available in full Stainless Steel / Stainless Steel top & Epoxy powder coated
- ACCESSORIES: Lithotomy bars, Stirrups, Step Stool & Mattress.

PATIENT TRANSFER EQUIPMENTS

21 (A) EMERGENCY & RECOVERY TROLLEY (CASUALTY TROLLEY)

Manually Height Adjustable :



Overall.Approx. Size :

L 1980mm x W 710 mm x Height adjustable from 700mm to 890mm

Features:

- Two sections removable & X-Ray permiable stretcher top.
- X-Ray cassette holder.
- Height adjustment & Trendelenburg / Reverse Trendelenburg position by screw & lever mechanism.
- Oxygen cylinder cage.
- I.V. Rod.
- S.S. swing away type side railings.
- Four heavy duty swivel castors of 125mm Dia. two with brakes
- Buffers at four corners.
- Pretreated & Epoxy powder coated.

21 (B) EMERGENCY & RECOVERY TROLLEY Height Adjustable by Hydraulic Pump



Overall.Approx. Size :

L 1980mm x W 710 mm x Height adjustable from 660mm to 910mm

Features:

- Two sections removable & X-Ray permiable stretcher top.
- X-Ray cassette holder.
- Height adjustment by foot operated Hydraulic pump.
- Gas spring assisted Trendelenburg / Reverse Trendelenburg positions.
- Oxygen cylinder cage.
- I.V. Rod.
- S.S. swing away type side railings.
- Four heavy duty swivel castors of 125mm Dia. two with brakes
- Buffers at four corners.
- Pretreated & Epoxy powder coated.

PATIENT TRANSFER EQUIPMENTS

19 (A) STRETCHER TROLLEY (FOLDING)

Overall approx. Size : L1830mm X W 570mm X H800mm

Folding approx. Size : L1000mm X W 570mm X H800mm

Features:

- Both side folding Stainless Steel top fixed on CRCA tube frame.
- Base Frame : 50 mm x 25 mm CRCA Tube.
- Mounted on 150 mm Dia heavy duty castors.
- Pretreated & Epoxy powder coated.
- Also available in Stainless Steel Top & Stainless Steel Base.



20 STRETCHER ON TROLLEY

Overall approx. Size :
L 2030 mm x W 610 mm x H 870 mm

Features:

- Base Frame : 50 mm x 25 mm CRCA Tube.
- Removable Stainless Steel top fixed on CRCA tube frame.
- Mounted on 150 mm Dia heavy duty castors.
- Pretreated & Epoxy powder coated.
- Also available in Stainless Steel.



HEALTH CHECK

42 (B) EXAMINATION COUCH WITH GAS SPRING ASSISTED BACKREST

Overall approx. Size : L 1830mm x W 575mm x H 850mm

Features:

- Frame : CRCA Rectangular Tubes, CRCA sheet cabinets.
- Gas spring assistance backrest.
- Three drawers & three cabinets.
- Sliding step stool.
- 100 mm thick reversible foam mattress.
- B.P. Apparatus stand.
- Pretreated & Epoxy powder coated.
- OPTIONAL: Lock facility for remaining Drawers & Cabinets.



42 (A) EXAMINATION COUCH

Overall approx. Size : L 1830mm x W 575mm x H 850mm

Features:

- Frame : CRCA Rectangular Tubes, CRCA sheet cabinets.
- Four cabinets - One with lock & one drawer with lock.
- Sliding step stool & Writing pad.
- 100 mm thick reversible foam mattress.
- Pretreated & Epoxy powder coated.
- OPTIONAL: Lock facility for remaining Cabinets.



MI - 43 GYANEC EXAMINATION COUCH

Overall approx. Size : L 1830mm x W 575mm x H 850mm

Features:

- Frame : CRCA Rectangular Tubes, CRCA sheet cabinets.
- Three cabinets one with lock & one drawer with lock.
- 100 mm thick reversible foam mattress with 'U' cut.
- Examination Trolley with three trays.
- Sliding step stool & writing pad.
- Pretreated & Epoxy powder coated.
- OPTIONAL: Lock facility for remaining Cabinets.



HEALTH CHECK

42 (C) STANDARD EXAMINATION COUCH



Overall approx. Size : L 1830mm x W 575mm x H 850mm

Features:

- Frame : CRCA Rectangular Tubes.
- CRCA sheet cabinet.
- One cabinet & one drawer with lock.
- Sliding step stool & Writing pad.
- 100 mm thick reversible foam mattress.
- Pretreated & Epoxy powder coated.

44 (A) GYANEC EXAMINATION TABLE

Overall approx. Size : L 1820mm x W 560mm x H 815mm

Features:

- Frame : Press bend CRCA sheet Top with backrest adjustable on ratchet.
- Perinial "U" cut at foot end.
- CRCA tubular legs.
- Step stool & Mattress.
- Lithotomy bars with Stirrups.
- Pretreated & Epoxy powder coated.



44 (B) PLAIN EXAMINATION TABLE



Features:

Overall approx. Size : L 1820mm x W 560mm x H 815mm

- Frame : Press bend CRCA sheet Top with backrest adjustable on ratchet.
- CRCA tubular legs.
- Step stool & Mattress.
- Pretreated & Epoxy powder coated.

O.T. EQUIPMENTS

28 INSTRUMENT TROLLEY

(A) INSTRUMENT TROLLEY Overall approx.
Size : L 710mm x B 460 mm x H 850 mm

(B) INSTRUMENT TROLLEY Overall approx.
Size : L 965mm x B 460 mm x H 850 mm

Features:

- Stainless Steel square tubes.
- Stainless Steel top & bottom shelves with 3 side railing on top shelf.
- Trolley mounted on 100 mm Dia. castors.
- Available in full Stainless Steel / powder coated frame with Stainless Steel shelves.
- Also available in round tube frame.
Supplied in KDC.



29 MAYO'S TROLLEY (SINGLE BAR)

Tray Size : L 555 mm x W 420mm
Height Adjustment From 845mm to 1300mm

Features:

- Frame: Stainless Steel Rectangular & Square Tubes.
- Trolley mounted on 50mm Dia. castors.
- Available in full Stainless Steel / Powder coated frame with S.S.Tray.



30 MAYO'S TROLLEY (DOUBLE BAR)

Tray Size : L 770mm x W 510mm
Height Adjustment From 845mm to 1300mm

Features:

- Frame : Stainless Steel Rectangular & Square Tubes.
- Trolley mounted on 50mm Dia. castors.
- Available in full Stainless Steel / Powder coated frame with S.S.Tray.

O.T. EQUIPMENTS



48 (A) WALL MOUNTED SCRUB

Size : W 482mm x D 457mm

Features:

- Length as per requirement .
- Fully Stainless Steel.
- Provided with instrument washing tray .
- OPTIONAL ACCESSORIES AT EXTRA COST:
 - a) Stainless Steel sleeper stand.
 - b) Long handle tap.

48 (B) FOOT OPERATED SCRUB

Overall approx .Size : W 558mm x D 457mm x H 1245mm

Features:

- Length as per requirement.
- Fully Stainless Steel.
- Instrument washing tray.



WARD REQUISITE



18 DRUG TROLLEY

A) Size : W 490mm x D 580mm x H 1040 mm,
4 trays & 100 mm Dia castors

B) Size : W 355mm x D 480mm x H 665mm

Features:

- 3 trays & 50mm Dia. castors.
- Pretreated & Epoxy powder coated.
- Available in Stainless Steel / Epoxy powder coated.
- Both are also available in Stainless Steel.

40 (A) BABY RESUSCITATION TROLLEY

Overall approx. Size : W 390mm x D 800mm x H 1080 mm

Features:

- Stainless Steel tray in head low position.
- Three drawer modular system.
- Adjustable I.V. Rod.
- Oxygen cylinder cage.
- Adjustable Examination Lamp.
- Bottom tray for mounting suction apparatus.
- Trolley mounted on 100 mm Dia castors.
- Pretreated & Epoxy powder coated.
- Available in Stainless Steel / Powder coated.



37 DRESSING TROLLEY WITH BOWL & BUCKET

Overall Trolley size : L 1080 x B 460 mm x H 850mm

Features:

- Frame : Stainless Steel square tubes of 25mm x 25mm.
- Two Stainless Steel shelves with railing on top.
- Stainless Steel bowl Dia 345 mm, bucket Dia 360mm.
- Trolley mounted on 100 mm Dia castors.
- Available in full Stainless Steel / Powder coated frame, with Stainless Steel shelves & railings.

WARD REQUISITE



31 REVOLVING STOOL

Overall approx. Size : Top 300mm Dia with
Height Adjustable 465 mm to 710 mm

Features:

- (A) Stainless Steel top & Epoxy powder coated base.
- (B) Stainless Steel Top & Stainless Steel base.
- (C) Stainless Steel Top & Stainless Steel base with Backrest.
- Optional : • Height adjustment by Gas spring • Cushions for seat & Backrest.



41 OXYGEN CYLINDER TROLLEY

Features:

- A) Pull Type with three castors 50mm Dia.
- B) Push Type with two castors 100mm Dia.
- C) Jumbo Oxygen Trolley with two castors with 125 mm Dia.

36 (A) TWO TIER BOWL STAND

Overall Approx Size : W 405mm x D 405 mm x H 840 mm

Features:

- Two Stainless Steel bowls Dia 360 mm each.
- Available in Stainless Steel.
- Epoxy powder coated.



34 SOILED LINEN TROLLEY

Linen Bag Size : 510mm dia x H 780mm

Features:

- Trolley mounted on 100mm Dia. three castors.
- Available in Stainless Steel/ Epoxy powder coated.



MI - 50 ECG MACHINE TROLLEY

Size: W355mm x D480mm x H665mm

Features:

- 3 trays, 50mm Dia. castors and Cable Hanger.
- Pre-treated & Epoxy powder coated.
- Available in Stainless Steel/Powder coated frame with S.S top.

MI - 35 KICK BUCKET

Overall Approx Size : W 365mm x D 365 mm x H 330 mm

Features:

- One stainless steel bowl Dia 360 mm.
- Available in Stainless Steel.
- Epoxy powder coated.



WARD REQUISITE



11 (A) CRIB ATTACHMENT



11 (B) CRIB WITH STAND



17 (C) FOOD TABLE-FULL S.S.
Stainless Steel frame and Stainless Steel top.
Overall approx Size : L 680 mm x W 360 mm x H 310 mm



10 (A) I.V. ROD FIX HEIGHT
Stainless Steel



10 (B) I.V. ROD TELESCOPIC
Stainless Steel



32 SALINE STAND
Ht. Adjustable from Min 1370 mm Max : 2435 mm
• Available in Full Stainless Steel / Epoxy powder coated stand & Stainless Steel saline Rod.



16 BED SIDE SCREEN

Overall approx. Size L open 2360mm /
closed 1200mm x W 550mm x H 1730mm

- CRCA tubular structure mounted on 50mm dia wheels two sides folding.



15 MULTIPURPOSE STOOL

Overall Top 280mm x 280mm & H 445mm
Frame : 25mm x 25mm, CRCA Square Tube

- Stainless Steel Top & Epoxy powder coated Base.
- Also Available in Stainless Steel Top & Stainless Steel base.



46 BED ELEVATOR

WHO-GMP, ISO 9001:2015, ISO 13485:2016 CERTIFIED COMPANY



Division of - Gajham Exports Pvt. Ltd.

Government Recognised Export House

"Ram House", 4, Gaiwadi Industrial Estate, S.V. Road, Goregaon (W),
Mumbai-400 062, INDIA. Tel.: +91-22-2879 4033 / 34 • Fax : +91-22-2879 4035
info@kinternational.co.in • www.kinternational.co.in